


# SURELOCK



## BONE MARROW BIOPSY SYSTEM

The SURELOCK needle has a sharp outer cannula that reduces the necessary penetration force. A turning motion instead of a twisting motion can be used to extract the sample in the inner stylet which significantly reduces patient trauma and insures a large intact bone marrow biopsy sample.

- Trapping system significantly reduces patient trauma
- Measuring stylet assures sample size before extraction
- Tapered tip traps specimen and allows for increased sampling size

Available sizes: 8 - 13G, 65 - 150mm

LEADERSHIP THROUGH IMAGINATION



**TSK LABORATORY, JAPAN**

製造販売 **株式会社 タスク**

〒328-0002 栃木県栃木市惣社町1485-11

TEL 0282-27-8426 FAX 0282-27-1943

ホームページ: [www.tsklab.co.jp](http://www.tsklab.co.jp) E-mail: [tsk@tsklab.co.jp](mailto:tsk@tsklab.co.jp)

Distributed by

